


### MIX. MOLD. DELIVER.

The CERAMENT® Bead Tray provides an adjunctive local delivery option for surgeons desiring targeted delivery to osseous defects and for dead space management.

- Made of flexible, medical-grade thermoplastic, making it easy to fill, set and release beads
- Three distinct bead sizes provide ultimate control over targeted delivery
- Supplied sterile and intended for single-use
- 10mL of CERAMENT will fill all bead cavities

BEAD SIZE	VOLUME
3mm	72 beads
4.7mm	60 beads
6mm	30 beads

CERAMENT Bead Tray is a delivery device for **CERAMENT® G** with Gentamicin, **CERAMENT® V** with Vancomycin, **CERAMENT® BONE VOID FILLER**



### MIX. MOLD. DELIVER.

#### Directions for use:

- Mix CERAMENT G, CERAMENT V, or CERAMENT BONE VOID FILLER according to their instructions for use
- Extrude the bone graft substitute paste into the bead mold
- Using the spatula, spread the material evenly, completely filling each bead cavity
- Remove air bubbles by gently tapping the bead mold on a flat surface
- Let the beads set for 15 minutes before releasing
- Flex the mold to release the beads
- Gently pack the beads into the treatment site
- Avoid overfilling or compressing the treatment site
- Close the wound(s) meticulously to avoid leakage into soft tissue



DESCRIPTION	ARTICLE NUMBER
CERAMENT® Bead Tray	A0513
CERAMENT® G 5mL	A0535-02
CERAMENT® G 10mL	A0535-01
CERAMENT® V 10mL	A0451-03
CERAMENT® BONE VOID FILLER 5mL	A0210-06
CERAMENT® BONE VOID FILLER 10mL	A0210-05
CERAMENT® BONE VOID FILLER 18mL	A0210-10