Čeština	3
Dansk	7
Deutsch	11
English	15
Ελληνικά	19
Español	23
Français	27
Hrvatski	31
Italiano	35
Nederlands	39
Norsk	43
Polski	47
Svenska	51
Українська	55

# CERAMENT™ | BONE VOID FILLER

CS

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## INSTRUCTIONS FOR USE

### DEVICE DESCRIPTION

CERAMENT™|BONE VOID FILLER is a fast-setting, injectable and moldable ceramic bone graft substitute intended for filling bone voids/gaps. The material consists of a powder and a liquid component. The major constituents of the powder are hydroxyapatite and calcium sulfate hemihydrate. The liquid component (C-TRU) contains Iohexol as a radio-opacification enhancer. Mixing the components, with the combined mixing injection (CMI) device, results in a viscous material intended to set *ex-vivo* or *in-vivo*. By combining hydroxyapatite and calcium sulfate an optimal balance is achieved between implant resorption rate and bone in-growth rate. Calcium sulfate acts as a resorbable carrier for hydroxyapatite. Hydroxyapatite has a slow resorption rate, high osteoconductivity promoting bone in-growth and gives long term structural support to newly formed bone.

The ceramic bone graft substitute is placed into the bone defect under visual inspection or under radiographic monitoring during open or percutaneous surgery.

The paste may be injected into the defect, molded by hand and digitally placed into the defect, or used to prepare beads that are placed into the defect. The accompanying injection device (ID) and Tip Extenders may be used to facilitate the filling of the bone defect.

When fully set *in vivo*, CERAMENT™|BONE VOID FILLER is drillable and can be used to augment hardware during the surgical procedure.

### PERFORMANCE

**Injectable:** 3-5 minutes.

**Molding by hand:** Initiation between 7-9 minutes.

**Beads:** Fill mold cavities before 5 minutes. Release beads after 15 minutes.

**Drillable:** After 15 minutes.

All times are from the start of mixing. More detailed information is provided in "Directions for use".

### INTENDED USE

CERAMENT™|BONE VOID FILLER is a resorbable ceramic bone graft substitute intended to fill gaps and voids in the skeletal system to promote bone healing.

CERAMENT™|BONE VOID FILLER provides a void/gap filler that during the surgical procedure can augment hardware and bone alignments.

### INDICATIONS

CERAMENT™|BONE VOID FILLER is indicated to be placed into bone voids or gaps in the skeletal system, i.e. extremities and pelvis (only during acetabular revision) not intrinsic to the stability of the bony structure. These osseous defects may be:

*In skeletally mature patients:* spontaneously occurring, surgically created, resulting from traumatic injury to the bone, identified during primary surgery and revision surgery, or osseous defects identified around hardware devices.

### CONTRAINDICATIONS

- Hypersensitivity to Iohexol or to any of the excipients
- Local infection at the site of implantation
- Pregnancy
- Breastfeeding
- Manifest thyroxicosis

### PRECAUTIONS

- Adhere to sterile surgical technique.

### Supportive therapy

- Control active bleeding and remove blood clots and tissue fragments if open surgery.
- Consult and comply with the IFU of any additional utensils.

### Patient related

- Pre existing calcium metabolism disorder (e.g. hypercalcemia).

### Device related

- Contact between CERAMENT™|BONE VOID FILLER and living bone, and provision of normal conditions of fracture healing or bone growth, are prerequisites for good treatment outcome.
- May cause inflammatory reaction if present in soft tissue.
- Not intended for load bearing areas, unless it can be assumed after thorough examination that the cortical bone, surgical fixation, or in situ osteosynthesis is sufficient for load bearing function.
- CERAMENT™|BONE VOID FILLER should not be used in joint spaces since its presence may cause inflammation or mechanical obstruction. When used in the vicinity of joints, care should be taken to ensure that bone voids implanted with CERAMENT™|BONE VOID FILLER are not in continuity with joint spaces, for example by ensuring adequate fracture reduction and/or by using fluoroscopic guidance (as determined by best surgical practice).
- Contact with joint fluid may cause resorption of CERAMENT™|BONE VOID FILLER.
- Overpressurization during injection should be avoided as intra-medullary injection with any bone void filler may lead to fat embolization or embolization of CERAMENT™|BONE VOID FILLER into the blood stream.
- Overpressurization of the device may lead to extrusion of the device beyond the site of its intended application and damage the surrounding tissues.
- Inject carefully under visual inspection or radiographic monitoring to avoid leakage.
- Do not overfill.
- In Aneurysmal bone cysts (ABCs) and other bone cysts prone to producing large volumes of fluid, there is increased risk of wound drainage, soft-tissue inflammation and wound breakdown if treated by open surgery. Use CERAMENT™|BONE VOID FILLER in bead form rather than complete void filling for these indications.
- No clinical experience with additives in CERAMENT™|BONE VOID FILLER. Using alternative mixing solutions and/or adding other substances to the mixture may affect the product setting in an uncontrolled manner, and may affect the safety and effectiveness of the product.
- When using as alternative to autograft, be aware that CERAMENT™|BONE VOID FILLER is not osteoinductive.
- If using CERAMENT™|BONE VOID FILLER in conjunction with allograft or autograft, apply each component separately, without intermixing before application, since intermixing may affect the setting in an uncontrolled manner.
- If using CERAMENT™|BONE VOID FILLER to augment hardware and bone alignments during surgical procedure, wait until material final setting time for optimal use.
- No clinical experience with prophylactic use.
- Do not use if the liquid is discolored or contains a precipitate.

### Related to general use of Iohexol

*These precautions have been established for the intravenous use of Iohexol-containing radiocontrast agents*

#### Hypersensitivity

- A positive history of allergy, asthma, or untoward reactions to iodinated contrast media indicates a need for special caution. Any application of contrast media should, therefore, be preceded by a detailed medical history, in patients with allergic diathesis and in patients with known hypersensitivity reactions a very strict indication is required.
- Premedication with corticosteroids or histamine H1 and H2 antagonists might be considered in patients at risk for intolerance, they may not, however, prevent anaphylactic shock, they may actually mask initial symptoms. In patients with bronchial asthma especially the risk for bronchospasm is increased.

en

- The risk of serious reactions in connection with use of Iohexol is regarded as minor. However, iodinated contrast media may provoke serious, life threatening, fatal anaphylactic reactions or other manifestations of hypersensitivity.
- Independent of quantity and route of administration, symptoms such as angio-oedema, conjunctivitis, coughing, pruritus, rhinitis, sneezing and urticaria may be indicative of a serious anaphylactoid reaction requiring treatment. A course of action should therefore be planned in advance, with necessary drugs and equipment available for immediate treatment, should a serious reaction occur. In imminent state of shock, administration of the contrast medium must be terminated immediately and - if necessary - specific intravenous treatment must be initiated.
- Patients using  $\beta$ -blockers may present with atypical symptoms of anaphylaxis which may be misinterpreted as vagal reaction.

#### Hydration

- Adequate hydration should be assured before and after contrast media administration. If necessary, the patient should be hydrated intravenously until excretion of the contrast medium is complete. This applies especially to patients with dys- and paraproteinaemias like multiple myeloma, diabetes mellitus, renal dysfunction, hyperuricaemia, as well as to elderly patients and patients in bad general condition. In patients at risk the water and electrolyte metabolism must be controlled and symptoms of a dropping serum calcium level must be taken care of. Due to the risk of dehydration induced by diuretics, at first, water and electrolyte rehydration is necessary to limit the risk of acute renal failure.

#### Cardio-circulatory reactions

- Care should also be taken in patients with serious cardiac disease/cardio-circulatory disease and pulmonary hypertension as they may develop hemodynamic changes or arrhythmias.
- Patients with cardiac insufficiency, severe coronary heart disease, instable angina pectoris, valvular diseases, previous myocardial infarction, coronary bypass and pulmonary hypertension are especially predisposed for cardiac reactions.
- In elderly patients and patients with pre-existing cardiac diseases, reactions with ischemic changes in the ECG and arrhythmia occur more frequently.

#### CNS disturbances

- Patients with acute cerebral pathology, tumors or a history of epilepsy are predisposed for seizures and merit particular care. Also alcoholics and drug addicts have an increased risk for seizures and neurological reactions.

#### Renal reactions

- Use of iodinated contrast media may cause contrast induced nephropathy, impairment of renal function or acute renal failure. To prevent these conditions following contrast media administration, special care should be exercised in patients with preexisting renal impairment and diabetes mellitus as they are at risk.
- Other predisposing factors are preceding renal failure following application of contrast media, a history of renal disease, age over 60 years, dehydration, advanced arteriosclerosis, decompensated cardiac insufficiency, high doses of contrast media and multiple injections, exposition to further nephrotoxins, severe and chronic hypertension, hyperuricaemia, paraproteinemias (myelomatosis and Waldenström's macroglobulinemia, plasmocytoma) or dysproteinemias.

#### Hepatic reactions

- A potential risk of transient hepatic dysfunction exists. Particular care is required in patients with severe disturbance of both renal and hepatic function as they may have significantly delayed contrast media clearance.

#### Myasthenia gravis

- The administration of iodinated contrast medium may aggravate the symptoms of myasthenia gravis.

#### Phaeochromocytoma

- In patients with phaeochromocytoma undergoing interventional procedures, alpha blockers should be given as prophylaxis to avoid hypertensive crisis.

#### Disturbed thyroid function

- Due to free iodide in the solutions and additional iodide released by deiodination, iodinated contrast media influence thyroid function. This may induce hyperthyroidism or even thyrotoxic crisis in predisposed patients.

- Patients with manifest but not yet diagnosed hyperthyroidism are at risk, patients with latent hyperthyroidism (e.g. nodular goitre) and patients with functional autonomy (often e.g. elderly patients, especially in regions with iodine deficiency) should therefore have their thyroid function assessed before examination if such conditions are suspected.
- Before administering an iodinated contrast agent, make sure that the patient is not about to undergo thyroid scan or thyroid function tests or treatment with radioactive iodine, as administration of iodinated contrast agents, regardless of the route, interferes with hormone assays and iodine uptake by the thyroid gland or metastases from thyroid cancer until urinary iodine excretion returns to normal.
- Following injection of an iodinated contrast agent, there is also a risk of induction of hypothyroidism.

#### Further risk factors

- Among patients with autoimmune diseases cases of serious vasculitis or Stevens-Johnson-like syndromes have been observed.
- Severe vascular and neurological diseases, especially in elderly patients are risk factors for reactions to contrast media.

#### SIDE EFFECTS

*The following complications have been reported to result from ceramic bone graft substitutes*

- Calcium based bone void fillers may color wound drainage white. This should not be a concern, however be aware of the risk of infection when drainage occurs.
- There have been reports in the literature on idiosyncratic reactions (laryngospasm and tachyarrhythmia) in children up to the age of 15 treated with ceramic bone substitute containing 75-100% calcium sulfate and 0-25% calcium phosphate.

#### INTERACTIONS

##### Related to general use of Iohexol

*These interactions have been established for the intravenous use of Iohexol-containing radioccontrast agents*

- There is a risk of the development of lactic acidosis when iodinated contrast agents are administered to diabetic patients treated with metformin, particularly in those with impaired renal function.
- Patients treated with Interleukin - 2 less than two weeks previously have been associated with an increased risk of delayed reactions (flu-like symptoms or skin reactions).
- The concomitant use of certain neuroleptics or tricyclic antidepressants can reduce the seizure threshold and thus increase the risk of contrast medium-induced seizures.
- Treatment with  $\beta$ -blockers may lower the threshold for hypersensitivity reactions, as well as necessitating higher doses of  $\beta$ -agonists when treating hypersensitivity reactions.
- Beta-blockers, vasoactive substances, angiotensin-converting enzyme inhibitors and angiotensin receptor antagonists may reduce efficacy of cardiovascular compensation mechanisms of blood pressure changes.
- All iodinated contrast media may interfere with tests on the thyroid function, thus the iodine binding capacity of thyroid may be reduced for up to several weeks.
- High concentrations of contrast media in serum and urine can interfere with laboratory tests for bilirubin, proteins or inorganic substances (e.g. iron, copper, calcium and phosphate).

#### DISCLAIMER

- In cases where it is not possible to establish a sufficient wound closure there might be a risk of skin inflammation reaction and/or prolonged wound drainage.
- Bone fracture and wound complications including hematoma, site drainage, infection and other complications are possible side effects of surgery.

## COMPONENTS, COMPOSITIONS & PARTS

**CERAMENT™ CMI**  
Mixing device pre-filled with ceramic bone graft substitute, a mixture of hydroxyapatite and calcium sulfate.

**CERAMENT™ C-TRU**  
Pre-filled syringe with iodine based mixing liquid. A water-soluble, radio-contrast enhancer (Iohexol) with the iodine concentration of 180 mg I/mL.

**Valve**  
To connect CERAMENT™ C-TRU and CERAMENT™ ID to CERAMENT™ CMI.

**CERAMENT™ ID**  
Injection device (accuracy of measuring scale  $\pm 5\%$ ). For the 18 mL product, two injection devices are included.

**Tip Extender**  
Two Tip Extenders in different length for use with CERAMENT™ ID facilitating paste injection.

## ADDITIONAL UTENSILS NEEDED

- Stopwatch

## OPTIONAL UTENSILS


- A cannula or needle with a minimum diameter of 16G
- Bead mold

## DIRECTIONS FOR USE

When handling CERAMENT™ BONE VOID FILLER adhere to sterile surgical techniques.

### Preparation of the paste - Step by step instructions

1. Retract the blue handle on the CERAMENT™ CMI and remove the red plunger stop ring.
2. Remove the plug on CERAMENT™ CMI and attach the Valve with the clear end to CERAMENT™ CMI by turning it clockwise.
3. Remove the plug from the liquid filled syringe (CERAMENT™ C-TRU) using the tool embedded in the blister tray bottom by pushing the syringe-plug into the tool and turn the syringe counter-clockwise. After removing the plug, press the CERAMENT™ C-TRU through the blue membrane of the Valve and attach it by turning it clockwise.
4. Empty the syringe with CERAMENT™ C-TRU into the CERAMENT™ CMI completely. Keep the plunger pushed to the bottom to avoid back flush. Detach CERAMENT™ C-TRU syringe from the Valve.
5. Start the stop watch ( $t = 0$  seconds). Mix in a tipped down position for 30 seconds with a frequency of approximately one complete stroke per second. Rotate the blue handle at the end positions.
6. Finish the mixing:
  - Fully retract the blue handle into its back position.
  - Lock the plunger by turning the blue collar clockwise 180° until a "click" is heard.
7. Attach the CERAMENT™ ID to the Valve by pressing it through the blue membrane and turning clockwise. Transfer the paste from the CERAMENT™ CMI immediately, with the grading on the CERAMENT™ ID facing the user. For the 18 mL product, the second syringe should be filled immediately after the first syringe. When the CERAMENT™ ID is completely filled, excess paste will begin to ooze from under the sleeve. Stop filling when this occurs.
8. Detach the filled CERAMENT™ ID from the Valve and remove the red plunger stop rod. The paste is now ready to use. If applicable, attach Tip Extender or an optional sterile needle (minimum 16G) to the CERAMENT™ ID.

## Filling of the bone void/gap - 3 different options

### Injection

- Wait until 3 minutes after start of mixing; carefully inject into the bone void/gap under visual inspection and/or by radiographic monitoring.
- Proceed until the void/gap is filled with an adequate amount of paste, as judged by the responsible physician.
- The paste is injectable between 3-5 minutes (from start of mixing) when using a 16 G needle.
- After the paste is *in situ*, allow to set for a few minutes before any adjustments are done or the wound is closed, especially if bleeding occurs.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

### Molded by hand

- Wait until approximately 3 minutes after start of mixing; extrude the paste to a solid mass of suitable size on sterile surface and wait until desired consistency for molding. The extrusion should be performed before 5 minutes.
- Molding of the paste by hand can be initiated between 7-9 minutes (the paste is moldable for a period of 1 minute).
- Place the molded product in the bone void or gap.
- Gently pack the molded material.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

### Molded with the use of Bead Mold Trays

- Select the size of beads needed.
- Fill the cavities in the mold before 5 minutes (from start of mixing).
- Allow paste to cure undisturbed for at least 15 minutes from start of mixing, before the beads are released from the mold.
- Flex mold to release the beads.
- Gently pack the beads in the bone void or gap.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

## STERILITY

CERAMENT™ BONE VOID FILLER is supplied sterile. The CERAMENT™ CMI is sterilized by gamma irradiation, the CERAMENT™ C-TRU is sterilized by steam and the surface sterilization of the complete device is by ethylene oxide.

The product is intended only for single use and cannot be re-sterilized by any method and shall not be re-used due to contamination risks.

## RESTRICTIONS

CERAMENT™ BONE VOID FILLER may only be sold and distributed for professional use.

Store CERAMENT™ BONE VOID FILLER unopened in a clean and dry environment at room temperature (15–30°C / 59–86°F).

The contents of this document may not be duplicated without written permission from BONESUPPORT AB.

Do not use if any of the packages are open or damaged or if the expiration date has been exceeded.

Excess material and opened but unused items must be discarded. Used material should be discarded in accordance with hospital procedures.

The products are protected by patents, see:  
<http://bonesupportpatents.com>

All trademarks herein are the property of BONESUPPORT AB, or its subsidiaries unless otherwise indicated.

Registered trademarks of BONESUPPORT AB:  
<https://www.bonesupport.com/trademarks.html>

## MANUFACTURED BY:

BONESUPPORT AB  
Scheelevägen 19  
IDEON Science Park  
SE-223 70 Lund, Sweden  
Email: info@bonesupport.com  
www.bonesupport.com

CE 2797

## EXPLANATION OF SYMBOLS

Do not use if package is damaged


Time, measured from start of mixing

## REF

A 0210-10 CERAMENT™ BONE VOID FILLER 18 mL  
A 0210-05 CERAMENT™ BONE VOID FILLER 10 mL  
A 0210-06 CERAMENT™ BONE VOID FILLER 5 mL


## Preparation of the paste


## Filling of the bone defect (void/gap) - 3 different options

All timelines are calculated from the Start of Mixing and vary slightly depending on product use


### Injection

### or Molded by hand\*

### or Molded with use of Bead Mold Trays\*\*

\* NOTE: Do not manipulate or touch extruded paste unnecessarily before initiation of molding.

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