Dansk	3
Deutsch	8
English	14
Español	19
Français	25
Italiano	31
Nederlands	37
Norsk	43
Polski	48
Suomi	54
Svenska	59

# CERAMENT™ |V

## BRUGSANVISNING

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## INSTRUCTIONS FOR USE

### DEVICE DESCRIPTION

CERAMENT™IV is a fast-setting, injectable and moldable ceramic bone graft substitute intended for the filling of bone voids/gaps. The material consists of powder and liquid components. The major constituents of the powder are hydroxyapatite, calcium sulfate hemihydrate and vancomycin. The liquid component (C-TRU) contains iohexol as a radio-opacification enhancer. Mixing the components, with the combined mixing injection device (CMI), results in a viscous material intended to set *ex-vivo* or *in-vivo*.

CERAMENT™IV delivers 66 mg vancomycin/mL paste.

By combining hydroxyapatite and calcium sulfate an optimal balance is achieved between implant resorption rate and bone in-growth rate. Calcium sulfate acts as a resorbable carrier for hydroxyapatite. Hydroxyapatite has a slow resorption rate, high osteoconductivity promoting bone in-growth and gives long term structural support to newly formed bone.

By adding vancomycin, colonization with vancomycin sensitive microorganisms can be prevented in order to protect bone healing.

The ceramic bone graft substitute is placed into the bone defect under visual inspection or under radiographic monitoring during open or percutaneous surgery.

The paste may be injected into the defect, molded by hand and digitally placed into the defect, or used to prepare beads that are placed into the defect. The accompanying injection device (ID) and Tip Extenders may be used to facilitate the filling of the bone defect.

When fully set *in vivo*, CERAMENT™IV can be used to augment hardware during the surgical procedure.

### PERFORMANCE

**Injectable:** approximately 3-5 minutes.

**Molding by hand:** Initiation between 8-12 minutes, and moldable for 1 minute.

**Beads:** Fill mold cavities before 5 minutes. Release beads after 20 minutes.

All times are from the start of mixing. More detailed information is provided in "Directions for use".

### INTENDED USE

CERAMENT™IV is a resorbable ceramic bone graft substitute intended to fill gaps and voids in the skeletal system to promote bone healing.

CERAMENT™IV provides a void/gap filler that during the surgical procedure can augment hardware and bone alignments.

Vancomycin is included in CERAMENT™IV to prevent colonization of vancomycin sensitive microorganisms in order to protect bone healing.

### INDICATIONS

CERAMENT™IV is indicated to be placed into bone voids or gaps in the skeletal system, i.e. extremities and pelvis (only during acetabular revision) not intrinsic to the stability of the bony structure. These osseous defects may be e.g. spontaneously occurring, surgically created, resulting from traumatic injury to the bone, identified during primary surgery and revision surgery, or osseous defects identified around hardware devices.

### CONTRAINDICATIONS

- Hypersensitivity to vancomycin
- Hypersensitivity to iohexol or any of the excipients (trometamol, sodium calcium edetate, hydrochloric acid), included in CERAMENT™(C-TRU)
- Pregnancy
- Breastfeeding
- Manifest thyroxicosis

### WARNINGS

Addition of CERAMENT™IV does not negate the need for systemic antibiotics.

### PRECAUTIONS

- Adhere to sterile surgical technique.
- Carefully examine patient medical history for relevant factors.

### Supportive therapy

- Consult and comply with the IFU of any additional used utensils.
- Control active bleeding and remove blood clots and tissue fragments if open surgery.
- Appropriate debridement is needed for successful use if pre-existing infection has been identified.
- Use of active suction drainage may lead to a decrease in effective dose of vancomycin.

### Device related

- Contact between CERAMENT™IV and living host bone, and provision of normal conditions of fracture healing or bone growth, are prerequisites for good treatment outcome.
- May cause inflammatory reaction if present in soft tissue.
- Not intended for load bearing areas, unless it can be assumed after thorough examination that the cortical bone, surgical fixation, or in situ osteosynthesis is sufficient for load bearing function.
- CERAMENT™IV should not be used in joints since its presence may cause inflammation or mechanical obstruction. When used in the vicinity of joints, care should be taken to ensure that bone voids implanted with CERAMENT™IV are not in continuity with joint spaces, for example by ensuring adequate fracture reduction and/or by using fluoroscopic guidance (as determined by best surgical practice).
- Contact with joint fluid may cause resorption of CERAMENT™IV.
- Overpressurization during injection should be avoided as intra-medullary injection with any bone void filler may lead to fat embolization or embolization of CERAMENT™IV into the blood stream.
- Overpressurization of the device may lead to extrusion of the device beyond the site of its intended application and damage the surrounding tissues.
- Inject carefully under visual inspection or radiographic monitoring to avoid leakage.
- Do not overfill.
- In Aneurysmal bone cysts (ABCs) and other bone cysts prone to producing large volumes of fluid, there is increased risk of wound drainage, soft-tissue inflammation and wound breakdown if treated by open surgery. Use CERAMENT™IV in bead form rather than complete void filling for these indications.
- No clinical experience with additives in CERAMENT™IV. Using alternative mixing solutions and/or adding other substances to the mixture may affect the product setting in an uncontrolled manner, and may affect the safety and effectiveness of the product.
- When used as an alternative to autograft, please note that like many synthetic bone void fillers CERAMENT™IV is not designed to be osteoinductive or osteogenic.
- If using CERAMENT™IV in conjunction with allograft or autograft, apply each component separately, without intermixing before application, since intermixing may affect the setting in an uncontrolled manner.
- If using CERAMENT™IV to augment hardware and bone alignments during surgical procedure, wait until material final setting time for optimal use.
- Do not use if the liquid is discolored or contains a precipitate.

### Patient related

- Pre existing calcium metabolism disorder (e.g. hypercalcemia).

### Related to general use of vancomycin

*These precautions have been established for the intravenous use of vancomycin*

- CERAMENT™IV should be used with caution in patients treated with systemic or topical tacrolimus, ciclosporin, muscle relaxants, loop diuretics and other potentially neurotoxic, nephrotoxic, or ototoxic drugs, such as amphotericin B, aminoglycosides, bacitracin, polymyxin B colistin, viomycin or cisplatin, and in patients with telcoplanin (glycopeptide antibiotic) sensitivity.
- Concomitant administration of vancomycin and anesthetic agents has been associated with erythema, histamine-like flushing and anaphylactoid reactions.

- Due to the potential ototoxicity and nephrotoxicity of vancomycin, CERAMENT™ IV should be used with care in patients with renal insufficiency. Regular monitoring of serum levels is advised in such patients.
- Concomitant use of vancomycin with muscle relaxants suxamethonium or vecuronium has been associated with the potential for enhanced neuromuscular blockade.
- In patients receiving vancomycin over a longer-term period or concurrently with other medications which may cause neutropenia or agranulocytosis, the leucocyte count should be monitored at regular intervals.
- Ototoxicity, which may be transitory or permanent has been reported in patients with prior deafness, who have received excessive intravenous doses, or who receive concomitant treatment with another ototoxic active substance such as an aminoglycoside. Vancomycin should also be avoided in patients with previous hearing loss.
- The elderly are particularly susceptible to auditory damage. Monitoring of vestibular and auditory function in the elderly should be carried out during and after treatment. Concurrent or sequential use of other ototoxic substances should be avoided.
- Stevens-Johnson syndrome (SJS) has been very rarely reported with the use of vancomycin. If symptoms or signs of SJS are present, specialized dermatological assessment should be sought.
- In case of severe persistent diarrhoea the possibility of pseudomembranous enterocolitis that might be life-threatening has to be taken into account. Anti-diarrhoeic medicinal products must not be given.

## Related to general use of iohexol

### *These precautions have been established for the intravenous use of iohexol-containing radiocontrast agents*

#### Hypersensitivity

- A positive history of allergy, asthma, or untoward reactions to iodinated contrast media indicates a need for special caution. Any application of contrast media should, therefore, be preceded by a detailed medical history, in patients with allergic diathesis and in patients with known hypersensitivity reactions a very strict indication is required.
- Premedication with corticosteroids or histamine H1 and H2 antagonists might be considered in patients at risk for intolerance, they may not, however, prevent anaphylactic shock, they may actually mask initial symptoms. In patients with bronchial asthma especially the risk for bronchospasm is increased.
- The risk of serious reactions in connection with use of iohexol is regarded as minor. However, iodinated contrast media may provoke serious, life threatening, fatal anaphylactic reactions or other manifestations of hypersensitivity.
- Independent of quantity and route of administration, symptoms such as angio-oedema, conjunctivitis, coughing, pruritus, rhinitis, sneezing and urticaria may be indicative of a serious anaphylactoid reaction requiring treatment. A course of action should therefore be planned in advance, with necessary drugs and equipment available for immediate treatment, should a serious reaction occur. In imminent state of shock, administration of the contrast medium must be terminated immediately and - if necessary - specific intravenous treatment must be initiated.
- Patients using  $\beta$ -blockers may present with atypical symptoms of anaphylaxis which may be misinterpreted as vagal reaction.

#### Hydration

- Adequate hydration should be assured before and after contrast media administration. If necessary, the patient should be hydrated intravenously until excretion of the contrast medium is complete. This applies especially to patients with dys- and paraproteinaemias like multiple myeloma, diabetes mellitus, renal dysfunction, hyperuricaemia, as well as to infants, small children, elderly patients and patients in bad general condition. In patients at risk the water and electrolyte metabolism must be controlled and symptoms of a dropping serum calcium level must be taken care of. Due to the risk of dehydration induced by diuretics, at first, water and electrolyte rehydration is necessary to limit the risk of acute renal failure.

#### Cardio-circulatory reactions

- Care should also be taken in patients with serious cardiac disease/cardiocirculatory disease and pulmonary hypertension as they may develop hemodynamic changes or arrhythmias.
- Patients with cardiac insufficiency, severe coronary heart disease, instable angina pectoris, valvular diseases, previous myocardial infarction, coronary bypass and pulmonary hypertension are especially predisposed for cardiac reactions.
- In elderly patients and patients with pre-existing cardiac diseases, reactions with ischemic changes in the ECG and arrhythmia occur more frequently.

#### CNS disturbances

- Patients with acute cerebral pathology, tumors or a history of epilepsy are predisposed for seizures and merit particular care. Also alcoholics and drug addicts have an increased risk for seizures and neurological reactions.

#### Renal reactions

- Use of iodinated contrast media may cause contrast induced nephropathy, impairment of renal function or acute renal failure. To prevent these conditions following contrast media administration, special care should be

exercised in patients with preexisting renal impairment and diabetes mellitus as they are at risk.

- Other predisposing factors are preceding renal failure following application of contrast media, a history of renal disease, age over 60 years, dehydration, advanced arteriosclerosis, decompensated cardiac insufficiency, high doses of contrast media and multiple injections, exposition to further nephrotoxins, severe and chronic hypertension, hyperuricaemia, paraproteinemias (myelomatosis and Waldenström's macroglobulinemia, plasmocytoma) or dysproteinemias.

#### Hepatic reactions

- A potential risk of transient hepatic dysfunction exists. Particular care is required in patients with severe disturbance of both renal and hepatic function as they may have significantly delayed contrast media clearance.

#### Myasthenia gravis

- The administration of iodinated contrast medium may aggravate the symptoms of myasthenia gravis.

#### Phaeochromocytoma

- In patients with phaeochromocytoma undergoing interventional procedures, alpha blockers should be given as prophylaxis to avoid hypertensive crisis.

#### Disturbed thyroid function

- Due to free iodide in the solutions and additional iodide released by deiodination, iodinated contrast media influence thyroid function. This may induce hyperthyroidism or even thyrotoxic crisis in predisposed patients.
- Patients with manifest but not yet diagnosed hyperthyroidism are at risk, patients with latent hyperthyroidism (e.g. nodular goitre) and patients with functional autonomy (often e.g. elderly patients, especially in regions with iodine deficiency) should therefore have their thyroid function assessed before examination if such conditions are suspected.
- Before administering an iodinated contrast agent, make sure that the patient is not about to undergo thyroid scan or thyroid function tests or treatment with radioactive iodine, as administration of iodinated contrast agents, regardless of the route, interferes with hormone assays and iodine uptake by the thyroid gland or metastases from thyroid cancer until urinary iodine excretion returns to normal.
- Following injection of an iodinated contrast agent, there is also a risk of induction of hypothyroidism.

#### Further risk factors

- Among patients with autoimmune diseases cases of serious vasculitis or Stevens-Johnson-like syndromes have been observed.
- Severe vascular and neurological diseases, especially in elderly patients are risk factors for reactions to contrast media.

#### Pediatric population

- Transient hyperthyroidism has been reported in premature infants, neonates and in other children after administration of iodinated contrast media. Premature infants are particularly sensitive to the effect of iodine. It is advisable to monitor thyroid function.
- Especially in infants and small children, adequate hydration should be assured before and after contrast media administration. Nephrotoxic medication should be suspended. The age dependent reduced glomerular filtration rate in infants can result in delayed excretion of contrast agent.
- Young infants (age < 1) and especially neonates are susceptible to electrolyte disturbance and hemodynamic alterations.

## SIDE EFFECTS

### *The following complications have been reported to result from ceramic bone substitutes*

- Calcium based bone void fillers may color wound drainage white. This should not be a concern, however be aware of the risk of infection when drainage occurs.
- There have been reports in the literature on idiosyncratic reactions (laryngospasm and tachyarrhythmia) in children up to the age of 15 treated with ceramic bone substitute containing 75-100% calcium sulfate and 0-25% calcium phosphate.

### *The following are complications related to the use of vancomycin*

- Ototoxicity and nephrotoxicity are the most common side effects associated with vancomycin therapy. Toxicity is related to serum concentration and it is advisable to check serum levels to confirm that the initial peak level does not exceed 26 micrograms/ml and that trough levels do not exceed 10 micrograms/ml.

*The following undesirable effects have been commonly associated with systemic administration of vancomycin and/or iohexol. The extent to which they apply to local orthopedic application is not known. They are listed here for information only.*

#### Vascular disorders

- Decrease in blood pressure

#### Respiratory, thoracic and mediastinal disorders

- Dyspnoea, stridor, transient changes in respiratory rate, respiratory distress

#### Skin and subcutaneous tissue disorders

- Flushing of the upper body ("red man syndrome"), exanthema and mucosal inflammation, pruritus, urticaria

### Renal and urinary disorders

- Renal insufficiency manifested primarily by increased serum creatinine and serum urea

### General disorders and administration site conditions

- Redness of the upper body and the face

*The following undesirable effects have been rarely, very rarely, or with unknown frequency associated with systemic administration of vancomycin and/or iohexol. The extent to which they apply to local orthopedic application is not known. They are listed here for information only.*

### Blood and the lymphatic system disorders

- Reversible neutropenia, agranulocytosis, eosinophilia, thrombocytopenia, pancytopenia

### Immune system disorders

- Hypersensitivity reactions (either immediate or delayed), anaphylactic reactions

### Ear and labyrinth disorders

- Transient or permanent loss of hearing, vertigo, tinnitus, dizziness

### Cardiac disorders

- Arrhythmia, bradycardia, myocardial infarction, cardiac failure, cyanosis, coronary artery spasm, chest pain, cardiac arrest

### Vascular disorders

- Hypertension, hypotension, vasculitis, flushing, shock, arterial spasm, thrombophlebitis, venous thrombosis

### Respiratory, thoracic and mediastinal disorders

- Cough, respiratory arrest, dyspnoea, severe respiratory symptoms and signs, pulmonary oedema, ARDS, bronchospasm, larygospasm, apnoea, aspiration, asthma attack

### Gastrointestinal disorders

- Nausea, pseudomembranous enterocolitis, vomiting, diarrhoea, acute pancreatitis or worsening of existing pancreatitis, salivary gland enlargement

### Skin and subcutaneous tissue disorders

- Exfoliative dermatitis, Stevens-Johnson syndrome, Lyell's syndrome, Linear IgA bullous dermatosis eosinophilia and systemic symptoms (DRESS syndrome), AGEP (Acute Generalized Exanthematous Pustulosis), rash, pruritus, urticaria, erythema multiforme, toxic epidermal necrolysis, acute generalised exanthematous pustulosis, drug rash with eosinophilia and systemic symptoms, psoriasis flare-up, erythema, drug eruption, skin exfoliation.

### Renal and urinary disorders

- Interstitial nephritis, acute renal failure, acute tubular necrosis, impaired renal function

### General disorders and administration site conditions

- Drug fever, shivering, pain and spasm of the chest and back muscles, pain and discomfort, administration site reactions, including extravasation, back pain asthenic conditions (including malaise, fatigue), feeling hot, hyperhidrosis, cold feeling, vasovagal reactions, pyrexia

### Injury, poisoning and procedural complications

- Iodism

### Nervous system disorders

- Headache, dysgeusia (transient metallic taste), syncope vasovagal, dizziness, paresis, paralysis, photophobia, somnolence, seizures, disturbance in consciousness cerebrovascular accident, sensory abnormalities (including hypoaesthesia), paraesthesia, tremor, transient motor dysfunction (including speech disorder, aphasia, dysarthria), transient contrast induced encephalopathy (including transient memory loss, coma, stupor, retrograde amnesia), disorientation, brain oedema

### Endocrine disorders

- Thyrototoxicosis, transient hypothyroidism

### Psychiatric disorders

- Confusion, agitation, restlessness, anxiety

### Eye disorders

- Visual impairment, Transient cortical blindness

## INTERACTIONS

### Related to general use of iohexol

*These interactions have been established for the intravenous use of iohexol-containing radiocontrast agents*

- There is a risk of the development of lactic acidosis when iodinated contrast agents are administered to diabetic patients treated with metformin, particularly in those with impaired renal function.
- Patients treated with Interleukin - 2 less than two weeks previously have been associated with an increased risk of delayed reactions (flu-like symptoms or skin reactions).

- The concomitant use of certain neuroleptics or tricyclic antidepressants can reduce the seizure threshold and thus increase the risk of contrast medium-induced seizures.

- Treatment with  $\beta$ -blockers may lower the threshold for hypersensitivity reactions, as well as necessitating higher doses of  $\beta$ -agonists when treating hypersensitivity reactions.

- Beta-blockers, vasoactive substances, angiotensin-converting enzyme inhibitors and angiotensin receptor antagonists may reduce efficacy of cardiovascular compensation mechanisms of blood pressure changes.

- All iodinated contrast media may interfere with tests on the thyroid function, thus the iodine binding capacity of thyroid may be reduced for up to several weeks.

- High concentrations of contrast media in serum and urine can interfere with laboratory tests for bilirubin, proteins or inorganic substances (e.g. iron, copper, calcium and phosphate).

## ADVERSE EVENTS

- If leakage occurs there is a risk of pulmonary embolization and in case of open foramen ovale (FO) additional risks in other organs.

## DISCLAIMER

- In cases where it is not possible to establish a sufficient wound closure there might be a risk of skin inflammation reaction and/or prolonged wound drainage.

- Bone fracture and wound complications including hematoma, site drainage, infection and other complications are possible side effects of surgery.

## COMPONENTS AND COMPOSITIONS

### CERAMENT™ CMI

Mixing device pre-filled with ceramic bone substitute, a mixture of hydroxyapatite and calcium sulfate.

### CERAMENT™ C-TRU

Pre-filled syringe with iodine based mixing liquid. A watersoluble, radio opacification enhancing component (iohexol) with the iodine concentration of 180 mg I/mL. To be used when preparing the vancomycin solution.

### Valve

To connect CERAMENT™ C-TRU and CERAMENT™ ID to CERAMENT™ CMI.

### CERAMENT™ ID

Injection device (accuracy of measuring scale  $\pm 5\%$ ).

### CERAMENT™ VANCOMYCIN

Glass vial with vancomycin. CERAMENT™ V will contain 66 mg vancomycin/mL paste.

### BONESUPPORT™ DP

A ventilated dispensing pin to facilitate easy handling when preparing the vancomycin solution.

### Tip Extender

Two Tip Extenders in different length to be used with CERAMENT™ ID facilitating paste injection.

## ADDITIONAL UTENSILS NEEDED

- Stopwatch

## OPTIONAL UTENSILS

- A cannula or needle with a minimum diameter of 16G
- Bead mold

## DIRECTIONS FOR USE

When handling CERAMENT™|V adhere to sterile surgical techniques.

### Dose recommendation

- 1 mL paste contains 66 mg vancomycin (provided as vancomycin hydrochloride).
- Adjust the volume of paste for appropriate dose of vancomycin to patient.
- Experience with local application of vancomycin: Topical usage of vancomycin dosages from 0.5 g up to 3 g in adult patients, without toxic serum levels or other safety issues have been described in peer reviewed publications.
- The recommendation for initial intravenous dose of vancomycin is:
  - *In patients aged 12 years and older:* The recommended dose is 15 to 20 mg/kg body weight every 8 to 12 hours (not to exceed 2 g per dose).
  - *In infants and children less than 12 years of age:* The recommended intravenous dose is 10 to 15 mg/kg body weight every 6 hours.
  - *In neonates and young infants:* the advice of a physician experienced in the management of neonates should be sought.


### Step by step instructions

#### Preparation of vancomycin solution

1. Remove the green cap from the CERAMENT™|VANCOMYCIN vial. Place the CERAMENT™|VANCOMYCIN vial on a solid surface and push the BONESUPPORT™|DP through the membrane of the vial.
2. Remove the plug from the liquid filled syringe (CERAMENT™|C-TRU) using the tool embedded in the blister tray bottom by pushing the syringe-plug into the tool and turn the syringe counter-clockwise. After removing the plug, attach the syringe to the BONESUPPORT™|DP.
3. Inject all CERAMENT™|C-TRU solution into the vial and dissolve the vancomycin powder.
4. Withdraw all of the vancomycin solution back into the CERAMENT™|C-TRU syringe. The vancomycin solution is now ready to be injected into the CERAMENT™|CMI, and should be used directly after having been prepared.

#### Preparation of the CERAMENT™|V paste

5. Retract the blue handle on the CERAMENT™|CMI and remove the red plunger stop ring.
6. Remove the plug of the CERAMENT™|CMI and attach the Valve with the clear end to CERAMENT™|CMI by turning it clockwise.
7. Attach the CERAMENT™|C-TRU syringe with vancomycin solution by pressing the syringe through the blue membrane of the Valve and turning the syringe clockwise.
8. Inject the vancomycin solution into the CERAMENT™|CMI. Empty the syringe completely. Keep the plunger pushed to the bottom to avoid back flush. Detach CERAMENT™|C-TRU syringe from the Valve.
9. Start the stop watch (t = 0 seconds). Mix in a tipped down position for 30 seconds with a frequency of approximately one complete stroke per second. Rotate the blue handle at the end positions.
10. Finish the mixing:
  - Fully retract the blue handle into its back position.
  - Lock the plunger by turning the blue collar clockwise until a "click" is heard.
11. Attach the CERAMENT™|ID to the Valve by pressing it through the blue membrane and turning clockwise. Transfer the paste immediately from the CERAMENT™|CMI with a light pressure, with the CERAMENT™|ID plunger pointing upwards and the grading facing the user. When the CERAMENT™|ID is completely filled, excess paste will begin to ooze from under the sleeve. Stop filling when this occurs.
12. Detach the filled CERAMENT™|ID from the Valve and remove the red plunger stop rod. The paste is now ready to use. If applicable, attach Tip Extender or an optional sterile needle (minimum 16G) to the CERAMENT™|ID.

## Filling of the bone void/gap - 3 different options

### Injection

- Wait until the paste has reached the desired consistency at approximately 3 minutes after start of mixing; carefully inject into the bone void/gap under visual inspection and/or by radiographic monitoring.
- Proceed until the void/gap is filled with an adequate amount of paste, as judged by the responsible physician.
- The paste is injectable between approximately 3-5 minutes (from start of mixing) when using a 16 G needle.
- After the paste is *in situ*, allow to set for a few minutes before any adjustments are done or the wound is closed, especially if bleeding occurs.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

### Molded by hand

- Wait until approximately 3 minutes after start of mixing; extrude the paste to a solid mass of suitable size on sterile surface and wait until desired consistency for molding. The extrusion should be performed before 5 minutes.
- Molding of the paste by hand can be initiated between 8-12 minutes (the paste is moldable for a period of 1 minute).
- Place the molded product in the bone void or gap.
- Gently pack the molded material.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

### Molded with the use of Bead Mold Trays

- Select the size of beads needed.
- Fill the cavities in the mold before 5 minutes (from start of mixing).
- Allow paste to cure undisturbed for at least 20 minutes from start of mixing, before the beads are released from the mold.
- Flex mold to release the beads.
- Gently pack the beads in the bone void or gap.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

## STERILITY

CERAMENT™|V is supplied sterile. The CERAMENT™|CMI is sterilized by gamma irradiation and the CERAMENT™|C-TRU is sterilized by steam. CERAMENT™|VANCOMYCIN is filtration sterilized and aseptically filled. Surface sterilization of the complete device is by ethylene oxide.

CERAMENT™|V is intended only for single use; the product cannot be re-sterilized by any method and shall not be re-used due to contamination risks.

## RESTRICTIONS

CERAMENT™|V may only be sold, distributed and used within the Intended use.

Store CERAMENT™|V unopened in a clean and dry environment at room temperature (15–25°C / 59–77°F).

The contents of this document may not be duplicated without written permission from BONESUPPORT AB.

Only for professional use.

Do not use if any of the packages are open or damaged or if the expiration date has been exceeded.

Excess material and opened but unused items must be discarded. Used material should be discarded in accordance with hospital procedures.

The products are protected by patents, see:

<http://bonesupportpatents.com>

All trademarks herein are the property of BONESUPPORT AB, or its subsidiaries unless otherwise indicated.

Registered trademarks of BONESUPPORT AB:  
<https://www.bonesupport.com/trademarks.html>

## MANUFACTURED BY:

BONESUPPORT AB  
Scheelevägen 19  
IDEON Science Park  
SE-223 70 Lund, Sweden  
Telephone: +46 46 286 53 70  
Fax: +46 46 286 53 71  
Email: [info@bonesupport.com](mailto:info@bonesupport.com)  
[www.bonesupport.com](http://www.bonesupport.com)

en

**REF**

A0451-03 CERAMENT™|V 10mL


**EXPLANATION OF SYMBOLS**

Do not use if package is damaged

Time, measured from start of mixing

en

**Preparation of the paste****Filling of the bone defect (void/gap) - 3 different options***All timelines are calculated from the Start of Mixing and vary slightly depending on product use***Injection****or Molded by hand\*****or Molded with use of Bead Mold Trays\*\***

\*NOTE: Do not manipulate or touch extruded paste unnecessarily before initiation of molding.

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Dokumentets innehåll får inte kopieras utan skriftligt medgivande från BONESUPPORT AB.
Alla varumärken i detta dokument tillhör BONESUPPORT AB eller dess dotterbolag, om inget annat anges.
Registrerade varumärken tillhörande BONESUPPORT AB: