MasterControl document

Document Number:	Revision:	Effective Date:
Title:		
CR/CO NO.	Previous Number:	(if from Paper original)
InfoCard Type:	Vault:	

Signature Manifest

Document Number: MC-001940

Revision: 01

Title: Multilingual IFU0027 for CERAMENT V A0581-01 CE

Effective Date: 27 Feb 2024

All dates and times are in GMT+1.

IFU0027 rev02 ML rev01 CV A0581-01 CE

Author

Name/Signature	Title	Date	Meaning/Reason
Carolina Wyrwa (CAROLINA.WYRWA)		13 Dec 2023, 03:18:07 PM	Approved

QA Approval

Name/Signature	Title	Date	Meaning/Reason
Therese Rick (THERESE.RICK)		13 Dec 2023, 03:26:45 PM	Approved

Quick Approval

Approve Now

Name/Signature	Title	Date	Meaning/Reason
Maria Bengzon (MARIA.BENZON)		27 Feb 2024, 03:29:09 PM	Approved

Quick Approval

Approve Now

Name/Signature	Title	Date	Meaning/Reason
Maria Bengzon (MARIA.BENZON)		27 Feb 2024, 03:29:20 PM	Approved

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CERAMENT™ V

INSTRUCTIONS FOR USE

DEVICE DESCRIPTION AND CLINICAL BENEFIT

CERAMENT V is an injectable and moldable ceramic bone graft substitute, consisting of calcium sulfate (60%), Hydroxyapatite (40%) and Vancomycin hydrochloride. The liquid component (CERAMENT C-TRU) contains iohexol (206 mg/mL paste, i.e. iodine concentration 95 mg/mL paste) as a radio-opacification enhancer. CERAMENT V delivers 66 mg Vancomycin/ mL paste. By combining Calcium sulfate and Hydroxyapatite an optimal balance is achieved between implant resorption rate and bone in-growth rate. Calcium sulfate acts as a resorbable carrier for Hydroxyapatite. Hydroxyapatite has a slow resorption rate, high osteoconductivity promoting bone in-growth and gives long term structural support to newly formed bone.

By adding Vancomycin, colonization with Vancomycin sensitive microorganisms can be prevented in order to protect bone healing.

INTENDED USERS

Medical professionals.

INTENDED PATIENT POPULATION

Skeletally mature patients who require surgery for bone voids which are at risk of bacterial colonization.

PERFORMANCE

The injecting performance is dependent upon the desired working consistency of the CERAMENT V.

Injectable: Approximately 3-5 minutes.

Molding by hand: Initiation between 8-12 minutes, and moldable for 1 minute.

Final setting: The wound can be closed at approximately 15 minutes. CERAMENT V attains final setting at approximately 25 minutes.

Beads: Fill mold cavities before 5 minutes. Release beads after 20 minutes.

Drillable: Drilling of the material can be performed at 15 minutes.

After implantation: CERAMENT V is resorbed and remodeled into new bone within 6-12 months after implantation.

More detailed information is provided in "Directions for use".

INTENDED USE

CERAMENT V is a resorbable ceramic bone graft substitute intended to fill gaps and voids in the skeletal system to promote bone healing. CERAMENT V provides a void/gap filler that during the surgical procedure can augment hardware and bone alignments.

Vancomycin is included in CERAMENT V to prevent colonization of Vancomycin sensitive microorganisms in order to protect bone healing.

INDICATIONS

CERAMENT V is indicated to be placed into bone voids or gaps in the skeletal system, i.e. extremities and pelvis (only during acetabular revision) not intrinsic to the stability of the bony structure. These osseous defects may be:

In skeletally mature patients: spontaneously occurring, surgically created, resulting from traumatic injury to the bone, identified during primary surgery and revision surgery, or osseous defects identified around hardware devices.

MODE OF ACTION

The device CERAMENT V has two modes of action:

- Primary mode of action is to be a resorbable ceramic bone graft substitute intended to fill gaps and voids in the skeletal system to promote bone healing. CERAMENT V provides a void/gap filler that during the surgical procedure can augment hardware and bone alignments.
- Secondary mode of action is to prevent colonization of Vancomycin sensitive microorganisms in order to protect bone healing.

CONTRAINDICATIONS

- Hypersensitivity to Vancomycin hydrochloride
- Hypersensitivity to iohexol or any of the excipients (Trometamol, Sodium calcium edetate, Hydrochloric acid), included in CERAMENT C-TRU
- Pregnancy
- Breastfeeding
- Manifest thyrotoxicosis

WARNINGS

The Vancomycin in CERAMENT V does not negate the need for systemic antibiotics.

PRECAUTIONS

Practical instruction for sterile / aseptic usage

- Adhere to sterile surgical technique when preparing the CERAMENT V paste.

Supportive therapy

- Consult and comply with the IFU of any additional used utensils.
- Control active bleeding and remove blood clots and tissue fragments if open surgery.
- Appropriate debridement is needed for successful use if pre-existing infection has been identified.
- Use of active suction drainage may lead to a decrease in effective dose of Vancomycin.

Device related

- Contact between CERAMENT V and vital bone, and provision of normal conditions of fracture healing or bone growth, are prerequisites for good treatment outcome.
- Not intended for load bearing areas, unless it can be assumed after thorough examination that the cortical bone, surgical fixation, or in situ osteosynthesis is sufficient for load bearing function.
- CERAMENT V should not be used in joints since its presence may cause inflammation or mechanical obstruction. When used in the vicinity of joints, care should be taken to ensure that bone voids implanted with CERAMENT V are not in continuity with joint spaces, for example by ensuring adequate fracture reduction and/or by using fluoroscopic guidance (as determined by best surgical practice).
- Contact with joint fluid may cause resorption of CERAMENT V.
- Overpressurization during injection should be avoided as intra-medullary injection with any bone void filler may lead to fat embolization or embolization of CERAMENT V into the blood stream.
- Do not overfill.
- In Aneurysmal bone cysts (ABCs) and other bone cysts prone to producing large volumes of fluid, there is increased risk of wound drainage, soft-tissue inflammation and wound breakdown if treated by open surgery. Use CERAMENT V in bead form rather than complete void filling for these indications.
- No clinical experience with additives in CERAMENT V. Using alternative mixing solutions and/or adding other substances to the mixture may affect the product setting in an uncontrolled manner, and may affect the safety and effectiveness of the product.
- When used as an alternative to autograft, please note that like any synthetic bone void fillers CERAMENT V is not designed to be osteoinductive or osteogenic.
- If using CERAMENT V in conjunction with allograft or autograft, apply each component separately, without intermingling before application, since intermingling may affect the setting in an uncontrolled manner.
- If using CERAMENT V to augment hardware and bone alignments during surgical procedure, wait until material final setting time for optimal use.
- Do not use if the liquid is discolored or contains a precipitate.

Patient related

- Pre existing calcium metabolism disorder (e.g. hypercalcemia).
- Careful examination of patient medical history is recommended.

Related to the use of Vancomycin

- CERAMENT V should be used with caution in patients treated with systemic or topical tacrolimus, ciclosporin, muscle relaxants, loop diuretics and other potentially neurotoxic, nephrotoxic, or ototoxic drugs, such as amphotericin B, aminoglycosides, bacitracin, polymyxin B colistin, viomycin or cisplatin, and in patients with teicoplanin (glycopeptide antibiotic) sensitivity.
- Concomitant administration of Vancomycin and anesthetic agents has been associated with erythema, histamine-like flushing and anaphylactoid reactions.
- Due to the potential ototoxicity and nephrotoxicity of Vancomycin, CERAMENT V should be used with care in patients with renal insufficiency. Regular monitoring of serum levels is advised in such patients.
- Concomitant use of Vancomycin with muscle relaxants suxamethonium or vecuronium has been associated with the potential for enhanced neuromuscular blockade.
- In patients receiving Vancomycin over a longer-term period or concurrently with other medications which may cause neutropenia or agranulocytosis, the leucocyte count should be monitored at regular intervals.
- Ototoxicity, which may be transitory or permanent has been reported in patients with prior deafness, who have received excessive intravenous doses, or who receive concomitant treatment with another ototoxic active substance such as an aminoglycoside. Vancomycin should also be avoided in patients with previous hearing loss.
- The elderly are particularly susceptible to auditory damage. Concurrent or sequential use of other ototoxic substances should be avoided.
- Stevens-Johnson syndrome (SJS) has been very rarely reported with the use of Vancomycin. If symptoms or signs of SJS are present, specialized dermatological assessment should be sought.

en

Related to the use of iohexol

Hypersensitivity

- A positive history of allergy, asthma, or untoward reactions to iodinated contrast media indicates a need for special caution. Any application of contrast media should, therefore, be preceded by a detailed medical history, in patients with allergic diathesis and in patients with known hypersensitivity reactions a very strict indication is required.
- Premedication with corticosteroids or histamine H1 and H2 antagonists might be considered in patients at risk for intolerance, they may not, however, prevent anaphylactic shock, which may actually mask initial symptoms. In patients with bronchial asthma especially the risk for bronchospasm is increased.
- The risk of serious reactions in connection with use of iohexol is regarded as minor. However, iodinated contrast media may provoke serious, life threatening, fatal anaphylactic reactions or other manifestations of hypersensitivity.
- Independent of quantity and route of administration, symptoms such as angio-oedema, conjunctivitis, coughing, pruritus, rhinitis, sneezing and urticaria may be indicative of a serious anaphylactoid reaction requiring treatment. A course of action should therefore be planned in advance, with necessary drugs and equipment available for immediate treatment, should a serious reaction occur. In imminent state of shock, administration of the contrast medium must be terminated immediately and - if necessary - specific intravenous treatment must be initiated.
- Patients using beta-adrenergic blocking agents, particularly asthmatic patients, may have a lower threshold for bronchospasm and are less responsive to treatment with beta agonists and adrenaline, which may necessitate the use of higher doses. These patients may also present with atypical symptoms of anaphylaxis which may be misinterpreted as vagal reaction.

Hydration

- Adequate hydration should be assured before and after contrast media administration. If necessary, the patient should be hydrated intravenously until excretion of the contrast medium is complete. This applies especially to patients with dys- and paraproteinaemias like multiple myeloma, diabetes mellitus, renal dysfunction, hyperuricaemia, as well as to elderly patients and patients in bad general condition. In patients at risk the water and electrolyte metabolism must be controlled and symptoms of a dropping serum calcium level must be taken care of. Due to the risk of dehydration induced by diuretics, at first, water and electrolyte rehydration is necessary to limit the risk of acute renal failure.

Cardio-circulatory reactions

- Care should also be taken in patients with serious cardiac disease/cardio-circulatory disease and pulmonary hypertension as they may develop hemodynamic changes or arrhythmias.
- Patients with cardiac insufficiency, severe coronary heart disease, instable angina pectoris, valvular diseases, previous myocardial infarction, coronary bypass and pulmonary hypertension are especially predisposed for cardiac reactions.
- In elderly patients and patients with pre-existing cardiac diseases, reactions with ischemic changes in the ECG and arrhythmia occur more frequently.

CNS disturbances

- Encephalopathy has been reported with the use of contrast media, such as iohexol. Contrast encephalopathy may manifest with symptoms and signs of neurological dysfunction such as headache, visual disturbance, cortical blindness, confusion, seizures, loss of coordination, hemiparesis, aphasia, unconsciousness, coma and cerebral oedema. Symptoms usually occur within minutes to hours after administration of iohexol, and generally resolve within days. Factors which increase blood-brain barrier permeability will ease the transfer of contrast media to brain tissue and may lead to possible CNS reactions for instance encephalopathy.

- Neurological symptoms caused by metastases, degenerative or inflammatory processes can be aggravated by application of contrast media.

Renal reactions

- Use of iodinated contrast media may cause contrast induced nephropathy, impairment of renal function or acute renal failure. To prevent these conditions following contrast media administration, special care should be exercised in patients with preexisting renal impairment and diabetes mellitus as they are at risk.
- Other predisposing factors are preceding renal failure following application of contrast media, a history of renal disease, age over 60 years, dehydration, advanced arteriosclerosis, decompensated cardiac insufficiency, high doses of contrast media and multiple injections, exposition to further nephrotoxins, severe and chronic hypertension, hyperuricaemia, paraproteinemias (myelomatosis and Waldenström's macroglobulinemia, plasmocytoma) or dysproteinemias.

Patients with disturbance of both hepatic and renal function

- Particular care is required in patients with severe disturbance of both renal and hepatic function as they may have significantly delayed contrast media clearance.

Myasthenia gravis

- The administration of iodinated contrast medium may aggravate the symptoms of myasthenia gravis.

Phaeochromocytoma

- In patients with phaeochromocytoma undergoing interventional procedures, alpha blockers should be given as prophylaxis to avoid hypertensive crisis.

Disturbed thyroid function

- Due to free iodine in the solutions and additional iodide released by deiodination, iodinated contrast media influence thyroid function. This may induce hyperthyroidism or even thyrotoxic crisis in predisposed patients.
- Before administering an iodinated contrast agent, make sure that the patient is not about to undergo thyroid scan or thyroid function tests or treatment with radioactive iodine, as administration of iodinated contrast agents, regardless of the route, interferes with hormone assays and iodine uptake by the thyroid gland or metastases from thyroid cancer until urinary iodine excretion returns to normal.
- Thyroid function tests indicative of hypothyroidism or transient thyroid suppression have been reported following iodinated contrast media administration.

Further risk factors

- Among patients with autoimmune diseases cases of serious vasculitis or Stevens-Johnson-like syndromes have been observed.
- Severe vascular and neurological diseases, especially in elderly patients are risk factors for reactions to contrast media.

Potential interactions with concomitant use of pharmaceutical agents, related to the use of Vancomycin

- Concurrent or sequential administration of Vancomycin with other potentially neurotoxic or/ and nephrotoxic active substances particularly gentamycin, amphotericin B, streptomycin, neomycin, kanamycin, amikacin, tobramycin, viomycin, bacitracin, polymyxin B, colistin, piperacillin/tazobactam and cisplatin may potentiate the nephrotoxicity and/ or ototoxicity of Vancomycin and consequently requires careful monitoring of the patient.
- Concurrent administration of Vancomycin and anaesthetic agents has been associated with erythema, histamine like flushing and anaphylactoid reactions.
- The effect (neuromuscular blockade) of muscle relaxants (such as succinylcholine) concurrently used can be enhanced and prolonged.

Potential interactions with concomitant use of pharmaceutical agents, related to the use of iohexol

- There is a risk of the development of lactic acidosis when iodinated contrast agents are administered to diabetic patients treated with metformin, particularly in those with impaired renal function.
- Patients treated with Interleukin - 2 less than two weeks previously have been associated with an increased risk of delayed reactions (erythema, flu-like symptoms or skin reactions).
- The concomitant use of certain neuroleptics or tricyclic antidepressants can reduce the seizure threshold and thus increase the risk of contrast medium-induced seizures.
- Treatment with β -blockers may lower the threshold for hypersensitivity reactions, as well as necessitating higher doses of β -agonists when treating hypersensitivity reactions.
- Beta-blockers, vasoactive substances, angiotensin-converting enzyme inhibitors and angiotensin receptor antagonists may reduce efficacy of cardiovascular compensation mechanisms of blood pressure changes.
- High concentrations of contrast media in serum and urine can interfere with laboratory tests for bilirubin, proteins or inorganic substances (e.g. iron, copper, calcium and phosphate).

SIDE EFFECTS

The following side effects have been reported to result from ceramic bone substitutes

- Calcium based bone void fillers may color wound drainage white. This should not be a concern, however be aware of the risk of infection when drainage occurs.
- May cause inflammatory reaction if present in soft tissue.
- There have been reports in the literature on idiosyncratic reactions (laryngospasm and tachyarrhythmia) in children up to the age of 15 treated with ceramic bone substitute containing 75-100% calcium sulfate and 0-25% calcium phosphate.

Known side effects related to the systemic use of Vancomycin

These side effects are reported related to systemic use of Vancomycin, but which are not known to be associated with the use of CERAMENT V. Nevertheless, it is recommended that users familiarize themselves with the potential side effects associated with the systemic use of Vancomycin.

- Ototoxicity and nephrotoxicity are side effects associated with systemic Vancomycin therapy; with toxicity related to serum concentration of Vancomycin.

Blood and the lymphatic system disorders

- Reversible neutropenia, agranulocytosis, eosinophilia, thrombocytopenia, pancytopenia.

Immune system disorders

- Hypersensitivity reactions, anaphylactic reactions

Ear and labyrinth disorders

- Transient or permanent loss of hearing, vertigo, tinnitus, dizziness

Vascular disorders

- Decrease in blood pressure, vasculitis

Respiratory, thoracic and mediastinal disorders

- Dyspnoea, stridor

Gastrointestinal disorders

- Nausea

Skin and subcutaneous tissue disorders

- Flushing of the upper body ("red man syndrome"), exanthema and mucosal inflammation, pruritus, urticaria

Renal and urinary disorders

- Renal insufficiency manifested primarily by increased serum creatinine and serum urea, interstitial nephritis, acute renal failure

General disorders and administration site conditions

- Phlebitis, redness of the upper body and the face, drug fever, shivering, pain and spasm of the chest and back muscles

Known side effects related to the systemic use of iohexol

These side effects are reported related to systemic use of iohexol, but which are not known to be associated with the use of CERAMENT V. Nevertheless, it is recommended that users familiarize themselves with the potential side effects associated with the systemic use of iohexol.

Immune system disorders

- Hypersensitivity reactions (either immediate or delayed), anaphylactic/anaphylactoid reactions

Nervous system disorders

- Headache

Gastrointestinal disorders

- Nausea, vomiting, abdominal pain

Cardiac disorders

- Bradycardia

General disorders and administration site conditions

- Feeling hot, hyperhidrosis, cold feeling, vasovagal reactions, pyrexia

DISCLAIMER

- In cases where it is not possible to establish a sufficient wound closure there might be a risk of skin inflammation reaction and/or prolonged wound drainage.
- Bone fracture and wound complications including hematoma, site drainage, infection and other complications are possible side effects of surgery.

COMPONENTS AND COMPOSITIONS

Preparation of Vancomycin solution

CERAMENT C-TRU

Pre-filled syringe with iodine based mixing liquid. A watersoluble, radio opacification enhancing component (iohexol) with the iodine concentration of 180 mg I/mL. To be used when preparing the Vancomycin solution.

CERAMENT VANCOMYCIN

Glass vial with Vancomycin hydrochloride. 1 mL of CERAMENT V will contain 66 mg Vancomycin.

BONESUPPORT DP

A ventilated dispensing pin to facilitate easy handling when preparing the Vancomycin solution.

Preparation of CERAMENT V paste

CERAMENT CMI

Mixing device pre-filled with ceramic bone substitute, a mixture of Calcium sulfate and Hydroxyapatite.

Valve

To enable connections between CERAMENT C-TRU and CERAMENT ID to CERAMENT CMI.

CERAMENT ID

Injection device (accuracy of measuring scale $\pm 5\%$).

Tip Extender

Two Tip Extenders in different lengths to be used with CERAMENT ID facilitating paste injection.

ADDITIONAL UTENSILS NEEDED

- Stopwatch

OPTIONAL UTENSILS

- A cannula or needle with a minimum diameter of 16G
- Bead mold

DIRECTIONS FOR USE

When handling CERAMENT V adhere to sterile surgical techniques.

Dose recommendation

- 1 mL paste contains 66 mg Vancomycin (provided as Vancomycin hydrochloride).
- Adjust the volume of paste for appropriate dose of Vancomycin to patient.
- Experience with local application of Vancomycin: Topical usage of Vancomycin dosages from 0.5 g up to 3 g in adult patients, without toxic serum levels or other safety issues have been described in peer reviewed publications.
- The recommendation for initial intravenous dose of Vancomycin is:
 - In patients aged 12 years and older: The recommended dose is 15 to 20 mg/kg body weight every 8 to 12 hours (not to exceed 2 g per dose).


Step by step instructions

CERAMENT V consists of an outer cardboard box containing a plastic tray in a Tyvek pouch (sterile barrier) containing all the components. The cardboard box with content is ethylene oxide sterilized to ensure surface sterility of all components.

Preparation of Vancomycin solution


- Remove the cap from the CERAMENT VANCOMYCIN vial. Place the CERAMENT VANCOMYCIN vial on a solid surface and push the BONESUPPORT DP through the membrane of the vial while keeping it stable. 
- Remove the plug from the liquid filled syringe (CERAMENT C-TRU) using the tool embedded in the blister tray bottom by pushing the syringe-plug into the tool and turn the syringe counter-clockwise. After removing the plug, attach the syringe to the BONESUPPORT DP. 
- Inject all CERAMENT C-TRU solution into the vial and dissolve the Vancomycin powder by gently swirling the vial. 
- Withdraw all of the Vancomycin solution back into the CERAMENT C-TRU syringe. The Vancomycin solution is now ready to be injected into CERAMENT CMI, and should be used directly after having been prepared. 

Preparation of the CERAMENT V paste

- Retract the blue handle on the CERAMENT CMI and remove the red plunger stop. 
- Remove the plug on CERAMENT CMI and attach the Valve with the clear end to the CERAMENT CMI by turning it clockwise. 
- Attach the CERAMENT C-TRU syringe with Vancomycin solution by pressing the syringe tip through the blue membrane of the Valve and turning the syringe clockwise. 
- Inject all Vancomycin solution into the CERAMENT CMI. Empty the syringe completely and avoid back flush by keeping the plunger pushed to the bottom before detaching it. 

en

- Detach the CERAMENT C-TRU syringe from the Valve on the CERAMENT CMI and immediately start the stop watch and start mixing (t=0 seconds). Mix in a tipped down position for 30 seconds with a frequency of approximately 1 complete stroke per second. Rotate the blue handle at the end positions.

- Fully retract the blue handle into its back position and lock the plunger by turning the blue collar clockwise until a "click" is heard.

- Attach the CERAMENT ID to the Valve and transfer the paste immediately with the CERAMENT ID plunger pointing upwards and with the grading on the CERAMENT ID facing the user. When the CERAMENT ID is completely filled, excess paste will begin to ooze from under the sleeve. Stop the transfer when this occurs.

- Detach the filled CERAMENT ID, remove the red plunger stopper and the paste is now ready to use. If applicable attach Tip Extender or an optional needle (minimum 16G) to the CERAMENT ID.

Filling of the bone void/gap - 3 different options

Injection

- Wait until the paste has reached the desired consistency at approximately 3 minutes after start of mixing; carefully inject into the bone void/gap under visual inspection and/or by radiographic monitoring.
- Proceed until the void/gap is filled with an adequate amount of paste, as judged by the responsible physician.
- The paste is injectable between approximately 3-5 minutes (from start of mixing) when using a 16 G needle.
- After the paste is *in situ*, allow to set for a few minutes before any adjustments are done or the wound is closed, especially if bleeding occurs.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

Molded by hand

- Wait until approximately 3 minutes after start of mixing; extrude the paste to a solid mass of suitable size on sterile surface and wait until desired consistency for molding. The extrusion should be performed before 5 minutes.
- Molding of the paste by hand can be initiated between 8-12 minutes (the paste is moldable for a period of 1 minute).
- Place the molded product in the bone void or gap.
- Gently pack the molded material.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

Molded with the use of Bead Mold Trays

- Select the size of beads needed.
- Fill the cavities in the mold before 5 minutes (from start of mixing).
- Allow paste to cure undisturbed for at least 20 minutes from start of mixing, before the beads are released from the mold.
- Flex mold to release the beads.
- Gently pack the beads in the bone void or gap.
- Close the wound(s) meticulously to avoid leakage into the soft tissue.
- Follow accepted clinical practice for postoperative care.

STERILITY

CERAMENT V is supplied sterile. The CERAMENT CMI is sterilized by gamma irradiation and the CERAMENT C-TRU is sterilized by steam. CERAMENT VANCOMYCIN is filtration sterilized and aseptically filled. Surface sterilization of the complete device is by ethylene oxide.

CERAMENT V is disposable and intended only for single use; the product should not be re-sterilized by any method and should not be re-used due to contamination risk.

INFORMATION TO BE GIVEN TO THE PATIENT

Implant card

The included implant card shall be provided to the patient. The information on the front side of the card should be completed and label A0550 be attached on the back side before it is handed over to the patient.

Patient leaflet

The information in this section shall be conveyed to the patient. The patient information is also available as a leaflet on the web page www.BONESUPPORT.com/patientinfo (indicated on implant card).

What is CERAMENT V?

CERAMENT V is a synthetic bone graft substitute, containing Calcium sulfate (60%), Hydroxyapatite (40%), Vancomycin hydrochloride and C-TRU. C-TRU is an iohexol containing liquid, which provides enhanced visibility on X-rays and under fluoroscopy. CERAMENT V contains iohexol (206 mg/mL paste), Vancomycin (66 mg/mL paste), trometamol (0.6 mg/mL paste), sodium calcium edetate (0.05 mg/mL paste) and hydrochloric acid (less than 1 mg/mL paste).

What is CERAMENT V used for?

CERAMENT V is used for patients that have a bone void/defect, in order to fill the bone void/defect to support bone healing and prevent infection.

How does CERAMENT V work?

Over time, CERAMENT V is resorbed and remodeled into new bone within 6-12 months after implantation. The Vancomycin helps to prevent colonization.

Is any special follow-up or monitoring required if you have CERAMENT V implanted?

No further surgical treatment or monitoring is required, you will have the routine follow-up appointments after your surgery.

Does the Vancomycin in CERAMENT V mean I don't have to take antibiotics separately?

No. If you have been prescribed antibiotics following your surgery, you must take them, as the Vancomycin in CERAMENT V does not replace the need for oral or intravenous antibiotics.

Adverse effects

Implanting CERAMENT V carries the same risks associated with every surgery such as infection, pain, bruising, swelling and bleeding at surgery site.

If CERAMENT V was used in proximity to a joint, this might cause an inflammation (redness, swelling, pain) of the joint.

Within soft tissue, CERAMENT V may cause an inflammatory reaction (redness, swelling, pain).

In rare cases CERAMENT V can lead to a collection of fluid or fluid that leaks from the wound, which is white in color. This "white wound drainage" resolves in most cases within two-three weeks without further treatment.

The risk of having a serious reaction to C-TRU, Vancomycin or their excipients (iohexol, trometamol, sodium calcium edetate, hydrochloric acid) is regarded as minor. However, swelling of the face, tongue, larynx, itchy or swollen skin, or a rash, can all be signs of an anaphylactic reaction.

Seek medical advice if you experience any of these symptoms.

Do you use other medicines?

Tell your surgeon if you are using, have recently used or might use any other medicines.

RESTRICTIONS

CERAMENT V may only be sold, distributed, and used in accordance with the intended use.

The contents of this document may not be duplicated without written permission from BONESUPPORT AB.

CERAMENT™ products are protected by patents:
<http://www.bonesupportpatents.com/>

CERAMENT™ is a registered trademark of BONESUPPORT AB.
<https://www.bonesupport.com/trademarks.html>

Store CERAMENT V unopened in a clean and dry environment at room temperature (15 – 25°C / 59 – 77°F).

Do not use if any of the packages are unintentionally opened or damaged or if the expiration date has been exceeded.

Excess material and opened but unused items must be discarded. Used material should be discarded in accordance with hospital procedures.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of safety and clinical performance (SSCP) will be available in the European database on medical devices (EUDAMED), where it is linked to the Basic UDI-DI number of the product, 0735005543046RZ.

<https://ec.europa.eu/tools/eudamed>

If you experience a serious incident in relation to CERAMENT V, report it immediately to **both** the manufacturer and the competent authority of your country.

MANUFACTURED BY:

BONESUPPORT AB
 Scheelevägen 19
 IDEON Science Park
 SE-223 70 Lund, Sweden
 T: +46 46 286 53 70
 Email: info@bonesupport.com
 Email: complaint@bonesupport.com (for incident reporting)
www.bonesupport.com

ORDERING INFORMATION

Art. Number / Item

REF

A0581-01 CERAMENT V 10 mL

CE 2797

en

Symbol glossary

Symbol	Description of symbol	Symbol	Description of symbol
	Time, measured from start of mixing		Do not use if package is damaged and consult IFU
	Unique Device Identifier (UDI)		Information website for patients
	Indicates that the device contains or incorporates a medicinal substance		Patient name or patient ID
	Indicates that the device is a medical device		Name and address of the implanting healthcare institution/ provider
	Indicates the sterile barrier system		Date of implantation

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Dokumentets innehåll får inte kopieras utan skriftligt medgivande från BONESUPPORT AB.
CERAMENT® är ett registrerat varumärke som tillhör BONESUPPORT AB: