

Below is a list of registered devices which can be used to provide an adjunctive delivery option for surgeons desiring percutaneous targeted delivery and controlled access to bone defects, subchondral defects, lesions and cysts.

DELIVERY CANNULAS

Drillable cannula provides complete command over placement while an ergonomic T-handle allows for easy control. Universal luer lock connector ensures easy fit for any size delivery syringe and the radiocontrast agent in CERAMENT BONE VOID FILLER provides good visibility for more precise graft placement.

Description	Part Number
BONESUPPORT Delivery Cannula** 11Ga x 120mm, closed tip, side port delivery	74389-01M
BONESUPPORT Delivery Cannula** 15Ga x 60mm, open tip, end port delivery	74388-01M

**Distributed by BONESUPPORT.

11Ga x 120mm,
Closed Tip,
Side Port Delivery



15Ga x 60mm,
Open Tip,
End Port Delivery



CERVOS CANNULAS

Drillable cannulas with streamlined design and attachable T-handle allows for more versatile delivery options.

Description	Part Number
CERVOS Access/Delivery 8Ga x 250mm**	CER-SUB-825
CERVOS Access/Delivery 11Ga x 110mm -- Open Tip**	CER-SUB-1111-OT

**Distributed by BONESUPPORT.

Drillable 11Ga x
110mm Open Tip

Drillable 8Ga x 250mm
Open Tip

Detachable T-handle



BONE MARROW HARVEST NEEDLE

Promote bone healing in cell-based orthopedic procedures. The double-diamond tip and side holes allow for efficient aspiration. An ergonomic twist-lock handle provides clinician comfort while the luer lock connector ensures a secure syringe attachment.

Description	Part Number
Bone Marrow Harvest Needle**	A0534-01

**Distributed by BONESUPPORT.



BEAD TRAY

Designed for targeted delivery of CERAMENT to osseous defects and for dead space management.

Description	Part Number
CERAMENT™ Bead Tray*	A0513

*Manufactured by BONESUPPORT.

Bead Volume & Size Guide	
3mm	72 beads
4.7mm	60 beads
6mm	30 beads



Manufacturer:

BONESUPPORT AB
Scheelevägen 19,
Ideon Science Park
SE- 223 70 Lund, Sweden
T: +46 46 286 53 70
F: +46 46 286 53 71
E: info@bonesupport.com

Ordering Contact:

T: 1.877.719.6718
E: us.sales@bonesupport.com

Complaint Contact:

E: complaint@bonesupport.com



BONESUPPORT, INC.,
60 William St, Suite 330
Wellesley, MA 02481

T: 1.877.719.6718
E: us.sales@bonesupport.com
W: www.bonesupport.com

