



THE SOLUTION FOR DELIVERING BIOLOGICS PRIOR TO FRACTURE FIXATION

- A 450mm cannula that allows for precise delivery of bone graft substitute within the medullary canal of long bones
- Side port cannula with a radiopaque tip allows for the precise delivery of biologics over a guide-wire after reaming
- Double cannulated design allows guide wire to remain in place during the delivery of biologics, without any concern of losing fracture reduction
- Holds ~11cc of bone substitute and goes over 3.2mm guide wire

PRESERVING LIMBS AND RESTORING HEALTH

2-CAN 450mm Delivery Cannula

2-CAN450B

