



THE SOLUTION FOR DELIVERING BIOLOGICS PRIOR TO FRACTURE FIXATION

- A 450mm cannula that allows for precise delivery of bone graft substitute within the medullary canal of long bones
- Side port cannula with a radiopaque tip allows for the precise delivery of biologics over a guide-wire after reaming
- Double cannulated design allows guide wire to remain in place during the delivery of biologics, without any concern of losing fracture reduction
- Holds 11.8cc of bone substitute
- 1Ga O.D. and 5Ga inner O.D. and goes over 3.2mm guide wire

PRESERVING LIMBS AND RESTORING HEALTH

2-CAN 450mm Delivery Cannula

2-CAN450B



Patented



Allows for the delivery of healing biologic materials and bone void filler while preventing problematic biologic extrusion.
www.Flow-FX.net

BONESUPPORT, INC.,
117 Fourth Ave, Suite 202
Needham, MA 02494

T: 1.877.719.6718
E: us.sales@bonesupport.com
W: www.bonesupport.com