



THE SOLUTION FOR DELIVERING BIOLOGICS PRIOR TO FRACTURE FIXATION

- A 350mm or 450mm cannula that allows for precise delivery of bone graft substitute within the medullary canal of long bones.
- Side port cannula with a radiopaque tip allows for the precise delivery of biologics over a guide-wire after reaming.
- Double cannulated design allows guide wire to remain in place during the delivery of biologics, without any concern of losing fracture reduction.



Patented

Product	Product Code	CC of CERAMENT® in Cannula	Gauges
2-CAN 450mm	2CAN450B	11.8cc	1Ga (8mm) outer and 5Ga inner (goes over 3.2mm or less K-wire)
2-CAN 350mm	2CAN350	9cc	1Ga (8mm) outer and 5Ga inner (goes over 3.2mm or less K-wire)

For complete product information, including indications, contraindications, warnings, precautions and potential adverse events, see package insert.